# **HOMEDICS®**

## **Body Manager**

### **Body Analyser Scale with USB PC Data Tracking Quick Start Guide**

---

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---

## GETTING STARTED

This is the Body Manager **Quick Start Guide**. To make full use of the body analyser features please download the full operating instructions from [www.body-manager.com](http://www.body-manager.com)

### **Step 1: Batteries**

Insert the 6 x AA batteries supplied into the battery compartment, ensuring the + and – terminals are correctly orientated.

### **Step 2: Carpet feet**

For use on carpet please attach the carpet feet provided (UK only).

### **Step 3: Software Installation**

Please go to [www.body-manager.com](http://www.body-manager.com) and download the latest version of the Body Manager software. The full operating instructions will be included with this download.

The Body Manager computer software is compatible with Windows 2000, XP and Vista. It is not compatible with Apple computers. To install the software follow this sequence:

1. Close all open programs.
2. Double-click on the application file that you have just downloaded.
3. Follow the on screen installation instructions.

### **Step 4: USB Flash Drive**

Before using the USB flash drive for the first time, insert it into the USB slot on your computer. Wait while the computer identifies the USB flash drive and installs the necessary driver.

To avoid corrupting the data on the USB flash drive always click the **Safely Remove Hardware** icon before removing the USB flash drive. The icon is on the right of the toolbar along the bottom edge of your computer screen. The icon is identified by a bright green arrow.

You are now ready to begin using your scale.

## SCALE OPERATION

This section explains how to take a simple weight reading. To make full use of the body analyser features please refer to the full operating instructions available from [www.body-manager.com](http://www.body-manager.com)

### **Weight Reading Only**

If you just want to know your weight, follow this sequence:

1. Press scale platform (on one of the corners nearest to you) firmly until you hear a beep then release.
2. Wait until display shows zero.
3. To change to st, kg or lb press button on underside of scale.
4. Stand on platform
5. Weight is displayed

**Note:** Weight only readings are not saved for transfer to your computer.

## INFORMATION FOR YOUR SAFETY

This product is for domestic use only.

**Professional medical guidance should always be sought before embarking on diet and exercise programmes.**

The Homedics Body Manager is calibrated for adult use (age 16+ years).

It is not suitable for use during pregnancy.

**CAUTION: Do not use this product if you have a pacemaker or other medical device fitted. If in doubt, consult your doctor.**

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**HoMedics Group Ltd**